




exoplan Supported Libraries

425 implant systems supported (7516 implants)



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Alpha-Bio Tec Ltd. (IL)							
ATID	13	●					● Ti, DC, PM
DFI	20	●					● Ti, DC, PM
I.C.E.	24	●					● Ti, DC, PM
NeO	24	●					● Ti, DC, PM
NeO-CS	14	●					● Ti, DC, PM
NICE	5	●					● Ti, DC, PM
SPI	24	●		FG not for Ø 6.0 mm			● Ti, DC, PM




FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available

MKT-exoplan-libraries-en, revision 2020-05-14



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Anthogyr SAS (FR)							
Axiom® BL 2.8	3						
Axiom® BL PX	21						● Ti, DC, PM
Axiom® BL REG	23						● Ti, DC, PM
Axiom® TL PX	83						
Axiom® TL REG	89						
 Argon Dental Vertriebs Gesellschaft mbH & Co. KG (DE)							
K3Pro® Rapid	28	●					● Ti, DC, PM
K3Pro® Short	12	●		FG not for Ø 6.0 mm			● Ti, DC, PM
K3Pro® Sure	20	●					● Ti, DC, PM

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 B&B Dental S.r.l. (IT)							
3P-SLIM	4						
3P-STANDARD	18						
EVOLUTION-SLIM	4						
EVOLUTION-STANDARD	15						
WIDE-STANDARD	20						
 BEGO Implant Systems GmbH & Co. KG (DE) 1/2							
Semados® Mini	9						
Semados® RI	19						● Ti
Semados® RS	26						● Ti
Semados® RSX	26						● Ti
Semados® S	26						● Ti
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 BEGO BEGO Implant Systems GmbH & Co. KG (DE) 2/2							
Semados® SC	27						● Ti
Semados® SCX	27						● Ti
 Bicon® , LLC (US)							
MAX 2.5™ Integra-CP™	6	●		FG not for Ø 4.0 mm, length 8 mm			● MU, Stock
Narrow® Integra-CP™	4	●		FG for Ø 3.0 mm			● MU, Stock
SHORT® Integra-CP™	10	●		FG not for Ø 6.0 mm			● MU, Stock
 BioHorizons® , Inc. (US) 1/2							
Internal RBT Surface with Laser-Lok®	16						
Internal RBT Surface	16						
Laser-Lok 3.0	3						
Mount-free Internal	16						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 BioHorizons® , Inc. (US) 2/2							
Mount-free Tapered Short	4						
Short Tapered Internal	4						
Tapered Immediate Molar	6						
Tapered Internal implants	24	●	●				
Tapered Plus	14	●	●				
Tapered Pro Implants	19	●	●				
Tapered Tissue Level	14	●	●				
 Bio-Micron (IT) 1/2							
BSC	5						
ICI SinOcta - engaging	10						● Ti, DC, PM, Stock
ICI SinOcta	14						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Bio-Micron (IT) 2/2							
IEI PEN	15						
ISC	8						
TWS_C	11						
TWS	13						
 Bionnovation® Produtos Biomédicos Ltda (BR) 1/2							
BioDIRECT Classic	20						
BioDIRECT Conic	22						
BioDIRECT SWE	25						
BioDIRECT XP	10						
BioINTERNO Conic	14						
BioINTERNO Medullary Conic	14						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Bionnovation® Produtos Biomédicos Ltda (BR) 2/2							
BioMORSE Conic	15						
BioMORSE EZ	15						
BioMORSE SWE	15						
BioMORSE XP	15						
Classic	25						
Conic	11						
 Biotech Dental (FR)							
KONTACT®	24						
WEEGO.fr®	14						



FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 bredent medical GmbH & Co. KG (DE)							
blueSKY	16						● Ti, DC, PM, MU
miniSKY	11						
narrowSKY	4						● Ti, DC, PM, MU
SKY classic	13						● Ti, DC, PM, MU
whiteSKY	13						● Ti
 BTI® Biotechnology Institute S.L. (ES) 1/2							
Externa Tiny®	31						
Externa® Universal	14						
Externa® Universal Plus	14						
Interna® 3.0	21						
Interna® Narrow	38						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 BTI® Biotechnology Institute S.L. (ES) 2/2							
Interna® Universal	37						
Interna® Universal Plus	23						
Interna® Wide	26						
 CAMLOG Biotechnologies GmbH (CH) 1/2							
CAMLOG® Progressive-Line Promote® Plus snap-in	15		●				● Ti, DC, PM, MU
CAMLOG® Progressive-Line Promote® Plus guide/screw-mounted	15						
CAMLOG® Root-Line 2 Promote® Plus	19		●				● Ti, DC, PM
CAMLOG® Screw-Line Promote® former-packaging	19		●				● Ti, DC, PM
CAMLOG® Screw-Line Promote® Plus Guided	11	●			● FG		● Ti, DC, PM
CAMLOG® Screw-Line Promote® Plus former-packaging	19		●				
CAMLOG® Screw-Line Promote® Plus snap-in	19		●				● Ti, DC, PM
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							




MKT-exoplan-libraries-en, revision 2020-05-14


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
camlog CAMLOG Biotechnologies GmbH (CH) 2/2							
CAMLOG® Screw-Line Promote® snap-in	19		●				● Ti, DC, PM
CERALOG®	6		●				
CONELOG® Progressive-Line Promote® Plus snap-in	19		●				● Ti, DC, PM
CONELOG® Progressive-Line Promote® Plus guide/screw-mounted	19						
CONELOG® Screw-Line Promote® Plus Guided	14	●			● FG		● Ti, DC, PM
CONELOG® Screw-Line Promote® Plus former-packaging	19		●				
CONELOG® Screw-Line Promote® Plus snap-in	19		●				● Ti, DC, PM
iSy®	12						● DC, PM
Conexão Conexão Sistemas de Prótese (BR) 1/2							
AR®	8						
Easy®	13						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Conexão Sistemas de Prótese Conexão Sistemas de Prótese (BR) 2/2							
Flash® Porous	10						
Flash® Vulcano®	10						
Grip® Porous	9						
Grip® Vulcano®	9						
Torq®	14						
 CONMET LLC (RU)							
Classic Implant	6						● Ti, DC
Conical Implant	15						● Ti, DC
Cylindrical Implant	25						● Ti, DC
Cylindrical Short Implant	3						● Ti, DC
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14




Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Cortex™ Dental Implants Industries Ltd. (IL)							
CLASSIX Conical	15						● Ti, DC, PM, MU
CLASSIX Internal Hex	15						● Ti, DC, PM
DYNAMIX Conical	19						● Ti, DC, PM, MU
DYNAMIX Internal Hex	15						● Ti, DC, PM
MAGIX Conical	15						● Ti, DC, PM, MU
 Cowellmedi Co., Ltd. (KR)							
INNO External	19						
INNO Internal Cuff 1.8 mm	19						
INNO Internal Cuff 2.4 mm	19						
INNO Narrow	8						
INNO Submerged	23						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 CSM Implant Co., Ltd. (KR)							
Submerged III	21						
Submerged	20						
 C-Tech Implant S.R.L. (IT)							
Esthetic Line	15	●		FG not for Ø 3.8 mm, length 8 mm			
Narrow Diameter	3	●					
 DENTAURUM GmbH & Co. KG (DE)							
tioLogic®	19						
tioLogic® ST	23						
tioLogic® TWINFIT	23						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Dentsply Sirona Inc. (US)							
Ankylos® C/X	22	●		FG not for Ø 3.5 mm + Ø 4.5 mm, lengths 6.6 mm + 17 mm; FG not for Ø 5.5 mm + Ø 7.0 mm	● FG		● Ti, DC, PM
Astra Tech OsseoSpeed™ EV	44						● Ti, DC, PM, MU
Astra Tech OsseoSpeed™ EV Profile	24						
Astra Tech OsseoSpeed™ TX	31						● Ti, DC, PM, MU
Astra Tech OsseoSpeed™ TX Profile	15						
Frialit®	17						● Ti, DC, PM
XiVE® S	25	●		FG not for Ø 3.4 mm, Ø 3.8 mm + Ø 4.5 mm, length 18 mm; FG not for Ø 5.5 mm	● FG		● Ti, DC, PM, MU
XiVE® TG	14						● DC, PM

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available


MKT-exoplan-libraries-en, revision 2020-05-14



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 etk (FR) <small>A LYRA COMPANY</small>							
Naturactis®	22						
Naturall+®	22						
 FGM Produtos Odontológicos (BR)							
Arcsys®	16						
 GC Corporation (JP)							
Aadva® Standard	11						
Aadva® Tapered	11						
GENESiO® Plus Straight	19						
SETiO® Plus Straight	19						




FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Implacil de Bortoli Material Odontológico Ltda (BR)							
CYLINDRICAL HE	20						
CYLINDRICAL HI	20						
TAPERED HE	15						
TAPERED HI	15						
TAPERED MORSE TAPER AR	20						
 Implant Direct™ Sybron International, LLC. (US) 1/3							
GoDirect®	15						● Ti
GoDirect Mini®	8						● Ti
InterActive™ - SBActive™	23						● Ti
InterActive™ - SBM™	23						● Ti
Legacy1™ - HA-coated	20						● Ti
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							





MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Implant Direct™ Sybron International, LLC. (US) 2/3							
Legacy1™ - SBM™	20						● Ti
Legacy2™ - HA-coated	40						● Ti
Legacy2™ - SBM™	40						● Ti
Legacy3™ - HA-coated	40						● Ti
Legacy3™ - SBM™	40						● Ti
Legacy4™ - HA-coated	40						● Ti
Legacy4™ - SBM™	40						● Ti
LegacyP™ - SBM™	31						● Ti
ScrewDirect® - SBM™	12						
ScrewIndirect® - SBM™	16						● Ti, MU
ScrewPlant® - HA-coated	15						● Ti
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							




Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Implant Direct™ Sybron International, LLC. (US) 3/3							
ScrewPlant® - SBM™	15						● Ti
ScrewPlus® - HA-coated	20						● Ti
ScrewPlus® - SBM™	20						● Ti
SwishActive™ - SBM™	16						● Ti
 Implant Logistics® Inc. (US)							
Implant One® 100 series	4						
Implant One® 300 series	12						
Implant One® 400 series	11						
Implant One® 500 series	10						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Keystone Dental, Inc. (US)							
Genesis	18						
Prima Plus™	15						
 KLOCKNER® Implant System (ES)							
ESSENTIAL® CONE	40	●		FG for Ø 3.5 mm, Ø 4.0 mm + Ø 4.5 mm; FG not for ESSENTIAL® CONE 1.5			● Ti, DC, PM
KL	48						● Ti, DC, PM
VEGA®	28						
VEGA®+	28						
 Leone® S.p.A. (IT) 1/2							
2.9	3						
Implant	12						
Max Stability	8						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



MKT-exoplan-libraries-en, revision 2020-05-14



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Leone® S.p.A. (IT) 2/2							
Short 6.5	1						
 LOGON OMF Bioactive Systems GmbH (DE)							
LOGON®	19						
 MaxiCare Implant System MaxiCare Dental Implants (SA)							
MaxiCare One Piece	8						
MaxiCare O-Ring	4						
MaxiCare Ortho	4						
MaxiCare V Implants	63						
 MEDENTIKA® A Straumann Group Brand MEDENTIKA® GmbH (DE) 1/2							
Microcone	30	●					● Ti, DC, PM
Minicone	3						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 MEDENTiKA® GmbH (DE) 2/2 <small>A Straumann Group Brand</small>							
Procone	11						● Ti, DC, PM
Quattrocone	19	●					● Ti, DC, PM
Quattrocone30	8						
 Medical Instinct® Deutschland GmbH (DE)							
BoneTrust® plus cone	19						● Ti, DC
BoneTrust® plus hex	19						
 MIS Implants Technologies Ltd. (IL) <small>Make It Simple</small>							
C1	19	●	● UNI	FG + PT for implant lengths 8 mm - 13 mm			● Ti, DC, PM
M4	30						● Ti, DC, PM, MU
SEVEN	30	●	● UNI	FG + PT for implant lengths 6 mm - 13 mm; FG not for Ø 6.0 mm			● Ti, DC, PM, MU
V3	19	●	● UNI	FG + PT for implant lengths 8 mm - 13 mm			● Ti
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ML Implant System® (AR)							
SHe®	19						
SHi®	22						
SHi Cortical®	15						
SRi®	15						
 NDI Medical Europe (DE) 1/2							
De Luxe	19						
Effective+	6						
Leader	19						
Optimum	19						
Prime	19						
Ridge+	2						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 NDI Medical Europe (DE) 2/2							
Short	3						
Standard Plus	15						
 Neobiotech Co., Ltd. (KR)							
IS-II active	19	●					
IS-III active	19	●					
IT-III active 1.8	19						
IT-III active 2.8	19						
IT-III active Wide Neck 1.8	15						
IT-III active Wide Neck 2.8	15						



FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
 DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Neodent® (BR) 1/2 <small>A Straumann Group Brand</small>							
Alvim CM Acqua®	15						
Alvim CM NeoPoros	15						
Drive® CM Acqua®	15						
Drive® CM NeoPoros	15						
Drive® GM® Acqua®	18						
Drive® GM® NeoPoros	18						
Facility Acqua®	3						
Facility NeoPoros	3						
Helix® GM® Acqua®	34						
Helix® GM® Long NeoPoros	6						
Helix® GM® NeoPoros	34						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Neodent® (BR) 2/2 <small>A Straumann Group Brand</small>							
Helix® HE™ Acqua®	30						
Titamax® CM NeoPoros	26						
Titamax® EX CM Acqua®	15						
Titamax® EX CM NeoPoros	17						
Titamax® GM® Acqua®	26						
Titamax® GM® NeoPoros	26						
Titamax® WS NeoPoros	6						



FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
 DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Neoss® GmbH (GB)							
ProActive® Straight Ø3.25	4	●		Only partially guided (guided drilling with manual implant insertion) available			
ProActive® Straight	24	●		FG and partially guided (guided drilling with manual implant insertion) available			
ProActive® Tapered	19	●		FG and partially guided (guided drilling with manual implant insertion) available			
ProActive® Wide	3	●		FG and partially guided (guided drilling with manual implant insertion) available			
 Nobel Biocare® Services AG (CH) 1/2							
Brånemark System® Mk III Groovy	22						
Brånemark System® Mk III Shorty	3						
Brånemark System® Mk III TiUnite®	24						
Brånemark System® Mk IV TiUnite®	13						
NobelActive® 3.0	4						
NobelActive®	24						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Nobel Biocare® Services AG (CH) 2/2							
NobelParallel™ Conical Connection	27						
NobelReplace® Conical Connection	15						
NobelReplace® Conical Connection PMC	15						
NobelReplace® Straight Groovy	14						
NobelReplace® Tapered Groovy	20						
NobelSpeedy® Groovy	37						
NobelSpeedy® Replace	19						
Replace Select™ Tapered PMC	20						
Replace Select™ Tapered TiUnite®	19						
Replace Select™ TC	8						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 OSSTEM IMPLANT CO., LTD (KR)							
SS Implant System	26						● DC, PM
TS Implant System	37	●		FG for implant lengths 7 mm - 13 mm (FG not for 7 mm extra short)			● Ti, DC, PM, MU
US Implant System	20						● Ti
 Paltop® Advanced Dental Solutions Ltd. (IL)							
Advanced Classic	21	●	●				● Ti, DC, PM, MU
Advanced +	21	●	●				● Ti, DC, PM, MU
DIVA	6	●	●				● Ti, DC, PM, MU
Dynamic	21	●	●				● Ti, DC, PM, MU
PAI - Paltop® Active Implant	21	●	●				● Ti, DC, PM, MU
PCA - Paltop® Conical Active	21	●	●				
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 SGS Dental Implant System Holding (CH)							
P1	27						
P1D	27						
P7	34						
P7D	34						
P7N	9						
 SIC invent AG (CH) 1/2							
SICace®	23						
SICmax® onepiece	2						
SICmax®	23						
SICtapered	20						
SICvantage® max	24						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 SIC invent AG (CH) 2/2							
SICvantage® tapered	24						
 S.I.N. - Sistema de Implante Nacional S.A. (BR) 1/2							
REVOLUTION®	4						● Ti
STRONG SW® CM	20						● Ti
STRONG SW® HE	23						● Ti
STRONG SW® HI	15						● Ti
STRONG SW® Plus CM	20						● Ti
STRONG SW® Plus HE	23						● Ti
STRONG SW® Plus HI	15						● Ti
TRYON®	23						● Ti
UNITITE® COMPACT	9						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 S.I.N. - Sistema de Implante Nacional S.A. (BR) 2/2							
UNITITE® PRIME	15						● Ti
UNITITE® SLIM	3						
 Southern Implants® (Pty) Ltd. (ZA) 1/2							
Deep Conical Cylindrical	18						● Ti, DC
Deep Conical Tapered	18						● Ti, DC
External Hex - Cylindrical	35		●	PT for implant lengths 8.5 mm - 15 mm; PT not for Ø 4.0 mm			
External Hex - Cylindrical-internal drive	26						
External Hex - MAX	11						
External Hex - MAX-8®	4						
External Hex - Platform Switch	5						
External Hex - Tapered	32	●	●				
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Southern Implants® (Pty) Ltd. (ZA) 2/2							
External Hex - Tapered-Internal Drive	14						
Internal Hex - Tapered	15						
IT-Cylindrical	19						
MAXIT®	11						
IT-Tapered	14						
PROMAX™	12						
Provata - Tapered	12	●		FG not for length 18 mm			
Tri-Nex - Cylindrical	15						● Ti, DC, PM, MU
TRI-MAX®	9						
Tri-Nex - Tapered	20						● Ti, DC, PM, MU
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
straumann group Straumann® Group (CH) 1/3							
Bone Level Roxolid® SLA® Loxim®	12						● Ti, DC, PM, Stock
Bone Level Roxolid® SLActive® Guided	9	●	●		● FG	● FG	● Ti, DC, PM, Stock
Bone Level Roxolid® SLActive® Loxim®	12						● Ti, DC, PM, Stock
Bone Level Tapered Roxolid® SLA® Loxim®	18						● Ti, DC, PM, Stock
Bone Level Tapered Roxolid® SLActive® Guided	9	●	●		● FG	● FG	● Ti, DC, PM, Stock
Bone Level Tapered Roxolid® SLActive® Loxim®	21						● Ti, DC, PM, Stock
Bone Level Tapered Titanium SLA® Loxim®	15						● Ti, DC, PM, Stock
Bone Level Titanium SLA® Loxim®	12						● Ti, DC, PM, Stock
Bone Level Titanium SLA®	12						● Ti, DC, PM, Stock
BLX - Bone Level X Roxolid® SLActive®	42						
Mini Implant	3						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
straumann group Straumann® Group (CH) 2/3							
PURE Ceramic Implant	4						● Ti
PURE Ceramic Monotype Implant	16						● Stock
Standard Plus Roxolid® SLA® Loxim®	22						● Ti, DC, PM, Stock
Standard Plus Roxolid® SLActive® Guided	9	●	●		● FG	● FG	● Ti, DC, PM, Stock
Standard Plus Roxolid® SLActive® Loxim®	22						● Ti, DC, PM, Stock
Standard Plus Short Roxolid® SLActive® Loxim®	3						● Ti, DC, PM, Stock
Standard Plus Titanium SLA® Loxim®	18						● Ti, DC, PM, Stock
Standard Plus Titanium SLA®	18						● Ti, DC, PM, Stock
Standard Roxolid® SLA® Loxim®	20						● Ti, DC, PM, Stock
Standard Roxolid® SLActive® Guided	6	●	●		● FG	● FG	● Ti, DC, PM, Stock
Standard Roxolid® SLActive® Loxim®	20						● Ti, DC, PM, Stock


FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available



MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
straumann group Straumann® Group (CH) 3/3							
Standard Titanium SLA® Loxim®	20						● Ti, DC, PM, Stock
Standard Titanium SLA®	20						● Ti, DC, PM, Stock
Tapered Effect Roxolid® SLA® Loxim®	11						● Ti, DC, PM, Stock
Tapered Effect Roxolid® SLActive® Guided	3	●	●		● FG	● FG	● Ti, DC, PM, Stock
Tapered Effect Roxolid® SLActive® Loxim®	11						● Ti, DC, PM, Stock
Tapered Effect Titanium SLA® Loxim®	11						● Ti, DC, PM, Stock
Tapered Effect Titanium SLA®	11						● Ti, DC, PM, Stock
sweden&martina Sweden & Martina S.p.A. (IT) 1/4							
CSR Reduced Neck	12						
CSR Straight Neck	19						
CSR Wide Neck	14						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							



MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Sweden & Martina S.p.A. (IT) 2/4							
Global	20						
Kohno One	8						
Kohno SP	11						
Kohno Straight	16						
Outlink ²	29						
Outlink ² SLIM	4						
Prama	18						
Prama machined	15						
Prama RF	17						
Prama RF machined	15						
Prama RF white	17						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							




Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Sweden & Martina S.p.A. (IT) 3/4							
Prama white	18						
Premium machined	22						
Premium Straight	24						
Premium Switching Platform machined	17						
Premium Switching Platform	17						
Premium TG	25						
Premium One	13						
Premium One machined	13						
Shelta	19						
Shelta SL	15						
Syra machined	15						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							





Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Sweden & Martina S.p.A. (IT) 4/4							
Syra Magnum	2						
Syra	15						
Syra Short	6						
Syra SL machined	16						
Syra SL	16						
 T.A.G. Dental (IL)							
Axis	25	●		FG for Ø 3.3 mm, Ø 3.75 mm + Ø 4.2 mm, lengths 8 mm - 13 mm			
Massif	21						


FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
 DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 TBR® Implants Group (FR)							
Baby 8	2						
Baby M	2						
Connect	15						
infinity	15						
M	15						
Z1-Connect	12						
Z1-infinity	15						
Z1-M	15						
 Thommen Medical AG (CH) 1/2							
CONTACT implant MC INICELL®	19	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 11.0 mm + 12.5 mm			● Ti, DC, PM, MU
CONTACT implant®	21	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 9.5 mm - 12.5 mm			● Ti, DC, PM, MU
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							


MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Thommen Medical AG (CH) 2/2							
CONTACT implant RC INICELL®	21	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 9.5 mm - 12.5 mm			● Ti, DC, PM, MU
ELEMENT implant LC INICELL®	11	●	●	FG for PF Ø 3.5 + PF Ø 4.0, lengths 6.5 mm - 12.5 mm			● Ti, DC, PM, MU
ELEMENT implant MC INICELL®	24	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 8.0 mm - 12.5 mm			● Ti, DC, PM, MU
ELEMENT implant®	30	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 6.5 mm - 12.5 mm			● Ti, DC, PM, MU
ELEMENT implant RC INICELL®	34	●	●	FG for PF Ø 3.5, PF Ø 4.0 + PF Ø 4.5, lengths 6.5 mm - 12.5 mm			● Ti, DC, PM, MU
 TRATE AG (CH)							
ROOTFORM	33						
 TRI® Dental Implants Int. AG (CH)							
TRI®-Narrow	4						● Ti, DC, PM, MU
TRI®-Octa	14						● Ti, DC, PM, MU
TRI®-Vent	17						● Ti, DC, PM, MU
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							


Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 Uris Inc. (US)							
OMNI System	28	●					
 VITA Zahnfabrik H. Rauter GmbH & Co. KG (DE)							
ceramic.implant CI	11						
<small>MaxiCare Implant System</small>  Vitane® Implant International (SA)							
Vitane®	63						
Vitane® One Piece	8						
Vitane® O-Ring	4						
Vitane® Ortho	4						
 Weihai WEGO Jericom Biomaterials Co., Ltd (CN)							
WEGO	16						● Ti, DC
<small>FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available</small>							

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ZIMMER BIOMET Zimmer Biomet (US) 1/5							
Certain® 3i T3® for Parallel Walled Navigator®	22						
Certain® 3i T3® with DCD® for Parallel Walled Navigator®	22						
Certain® Full OSSEOTITE® for Parallel Walled Navigator®	20						
Certain® NanoTite® for Parallel Walled Navigator®	22						
Certain® OSSEOTITE® for Parallel Walled Navigator®	23						
Certain® OSSEOTITE®2 for Parallel Walled Navigator®	23						
Certain® PREVAIL® 3i T3® for Parallel Walled Navigator®	15						
Certain® PREVAIL® 3i T3® with DCD® for Parallel Walled Navigator®	15						
Certain® PREVAIL® NanoTite® for Parallel Walled Navigator®	10						
Certain® PREVAIL® OSSEOTITE® for Parallel Walled Navigator®	10						
Certain® PREVAIL® OSSEOTITE®2 for Parallel Walled Navigator®	15						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							


MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ZIMMER BIOMET Zimmer Biomet (US) 2/5							
Certain® Tapered 3i T3® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered 3i T3® with DCD® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered CP4 TAPERED for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered Full OSSEOTITE® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered NanoTite® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered OSSEOTITE® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm			● Ti, DC, PM
Certain® Tapered PREVAIL® 3i T3® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm			
Certain® Tapered PREVAIL® 3i T3® with DCD® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm			
Certain® Tapered PREVAIL® CP4 TAPERED for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm			
Certain® Tapered PREVAIL® NanoTite® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm			
External Hex 3i T3®	34						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							


MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ZIMMER BIOMET Zimmer Biomet (US) 3/5							
External Hex 3i T3® with DCD®	30						
External Hex Full OSSEOTITE®	29						
External Hex NanoTite®	27						
External Hex OSSEOTITE®	34						
External Hex OSSEOTITE®2	30						
External Hex Tapered 3i T3®	20						
External Hex Tapered 3i T3® with DCD®	20						
External Hex Tapered Full OSSEOTITE®	20						
External Hex Tapered NanoTite®	20						
External Hex Tapered OSSEOTITE®	20						
Eztetic® - MTX® Surface, 0.5mm Machined Collar with Microgrooves	5						
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							













MKT-exoplan-libraries-en, revision 2020-05-14











Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ZIMMER BIOMET Zimmer Biomet (US) 4/5							
Eztetic® - MTX® Surface, Fully Textured with Microgrooves	5						
Spline® Twist - MP-1® HA Surface	11						
Spline® Twist - MTX® Surface	11						
Tapered Screw-Vent® - MP-1® HA Dual Transition Selective Surface	20						● Ti, DC, PM
Tapered Screw-Vent® - 0.5 mm Machined Collar, Microgrooves, MP-1® HA Dual Transition Selective Surface	20						● Ti, DC, PM
Tapered Screw-Vent® - 0.5 mm Machined Collar, MTX® Surface, Microgrooves	20						● Ti, DC, PM
Tapered Screw-Vent® - Full MTX® Surface Texturing, Microgrooves	20						● Ti, DC, PM
Tapered Screw-Vent® - MTX® Surface	20						● Ti, DC, PM
Tapered Screw-Vent® - MTX® Textured Collar, Microgrooves, MP-1® HA Dual Transition Selective Surface	20						● Ti, DC, PM
Tapered SwissPlus®	12						
Trabecular Metal™ - MTX® Surface, 0.5 mm Machined Collar with Microgrooves	13						● Ti, DC, PM
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available							

MKT-exoplan-libraries-en, revision 2020-05-14

Available Implant Systems		Available Original Sleeves			Compatible Sleeves		Prosthetics
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	Compatible prosthetic components
 ZIMMER BIOMET Zimmer Biomet (US) 5/5							
Trabecular Metal™ - MTX® Surface, Fully Textured with Microgrooves	13						● Ti, DC, PM

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available
 DC = Direct-to-Implant Interfaces available, Ti = Titanium Bases available, Stock = Stock Abutments available, MU = Multi-Unit Abutments available, PM = Pre-milled blanks available

Sleeve / Anchor Pin Manufacturer	Fully guided Sleeves Compatible with	Universal Sleeves Compatible with	Pilot Sleeves Compatible with	Anchor / Fixation Pins	
 3D Diagnostix •COM	3D Diagnostix, Inc. (US)		all implants	all implants	●
	Argon Dental Vertriebs Gesellschaft mbH & Co. KG (DE)	Argon Dental Vertriebs Gesellschaft mbH & Co. KG			
	Bicon®, LLC (US)	Bicon®, LLC			
	BioHorizons®, Inc. (US)	BioHorizons®, Inc.		BioHorizons®, Inc.	●
	CAMLOG Biotechnologies AG (CH)	CAMLOG Biotechnologies GmbH		CAMLOG Biotechnologies GmbH	
	C-Tech Implant S.R.L. (IT)	C-Tech Implant S.R.L.			●
	Dentsply Sirona Inc. (US)	Dentsply Sirona Inc.			
	KLOCKNER® Implant System (ES)	KLOCKNER® Implant System			●
	MEDENTiKA® GmbH (DE)	MEDENTiKA® GmbH			
	MIS Implants Technologies Ltd. (IL)	MIS Implants Technologies Ltd.		MIS Implants Technologies Ltd.	●
	Neobiotech Co., Ltd. (KR)	Neobiotech Co., Ltd.			●
	Neoss® GmbH (GB)	Neoss® GmbH (inclusive partially guided)			

Sleeve / Anchor Pin Manufacturer	Fully guided Sleeves Compatible with	Universal Sleeves Compatible with	Pilot Sleeves Compatible with	Anchor / Fixation Pins
 NT digital implant technology	nt-trading GmbH & Co. KG (DE)	CAMLOG Biotechnologies GmbH, Dentsply Sirona Inc., Straumann® Group		all implants
	OSSTEM IMPLANT CO., LTD (KR)	OSSTEM IMPLANT CO., LTD		●
	Paltop® Advanced Dental Solutions Ltd. (IL)	Paltop® Advanced Dental Solutions Ltd.		●
	Southern Implants® (Pty) Ltd. (ZA)	Southern Implants® (Pty) Ltd.		Southern Implants® (Pty) Ltd.
	steco-system-technik GmbH & Co. KG (DE)	Straumann® Group	all implants	all implants
	Straumann® Group (CH)	Straumann® Group		Straumann® Group
	T.A.G. Dental (IL)	T.A.G. Dental		●
	Thommen Medical AG (CH)	Thommen Medical AG		Thommen Medical AG
	Uris Inc. (US)	Uris Inc.		●
	Zimmer Biomet (US)	Zimmer Biomet		